

Rat Opa1 / Dynamin-like 120 kDa protein, mitochondrial ELISA Kit

Catalog No: #EK21084



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Rat Opa1 / Dynamin-like 120 kDa protein, mitochondrial ELISA Kit |
| Brief Description | ELISA Kit |
| Applications | ELISA |
| Species Reactivity | Rt |
| Other Names | Optic atrophy protein 1 homolog |
| Accession No. | Uniprot:Q2TA68 |
| Storage | 6 months |

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.054ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.