

Rat Slc1a2 / Excitatory amino acid transporter 2 ELISA Kit

Catalog No: #EK21085



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Slc1a2 / Excitatory amino acid transporter 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Other Names	GLT-1,Sodium-dependent glutamate/aspartate transporter 2,GLUT-R,Solute carrier family 1 member 2,Eaat2,Glt1
Accession No.	Uniprot:P31596
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.147ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.