

# Human PEX1 / Peroxisome biogenesis factor 1 ELISA Kit

Catalog No: #EK21086



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human PEX1 / Peroxisome biogenesis factor 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Peroxin-1, Peroxisome biogenesis disorder protein 1
Accession No.	Uniprot:O43933
Storage	6 months

## Application Details

Detection Range: 0.312-20ng/mL Assay Type: Sandwich Sensitivity: 0.135ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.