

Human PAWR / PRKC apoptosis WT1 regulator protein ELISA Kit

Catalog No: #EK21087



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human PAWR / PRKC apoptosis WT1 regulator protein ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Prostate apoptosis response 4 protein,Par-4,PAR4
Accession No.	Uniprot:Q96IZ0
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:35pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.