

Human P3H4 / Endoplasmic reticulum protein SC65 ELISA Kit

Catalog No: #EK21088



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human P3H4 / Endoplasmic reticulum protein SC65 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Leprecan-like protein 4, Nucleolar autoantigen No55, Prolyl 3-hydroxylase family member 4, LEPREL4, NOL55, SC65, Synaptonemal complex protein SC65
Accession No.	Uniprot:Q92791
Storage	6 months

Application Details

Detection Range: 78-5000pg/mL Assay Type: Sandwich Sensitivity: 35pg/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.