

# Human RPGR / X-linked retinitis pigmentosa GTPase regulator ELISA Kit

Catalog No: #EK21090



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human RPGR / X-linked retinitis pigmentosa GTPase regulator ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	RP3,XLRP3
Accession No.	Uniprot:Q92834
Storage	6 months

## Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.33ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.