

Human EXT1 / Exostosin-1 ELISA Kit

Catalog No: #EK21091



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human EXT1 / Exostosin-1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Glucuronosyl-N-acetylglucosaminyl-proteoglycan/N-acetylglucosaminyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase, Multiple exostoses protein 1, Putative tumor suppressor protein EXT1
Accession No.	Uniprot:Q16394
Storage	6 months

Application Details

Detection Range: 0.78-50ng/mL Assay Type: Sandwich Sensitivity: 0.38ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.