

Human ELF1 / ETS-related transcription factor Elf-1 ELISA Kit

Catalog No: #EK21093



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ELF1 / ETS-related transcription factor Elf-1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	E74-like factor 1
Accession No.	Uniprot:P32519
Storage	6 months

Application Details

Detection Range:0.312-20ng/ml Assay Type:Sandwich Sensitivity:0.095ng/ml

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.