

# Human EPRS / Bifunctional glutamate/proline--tRNA ligase ELISA Kit

Catalog No: #EK21094



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human EPRS / Bifunctional glutamate/proline--tRNA ligase ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Bifunctional aminoacyl-tRNA synthetase, Cell proliferation-inducing gene 32 protein, Glutamaryl-prolyl-tRNA synthetase, GLNS, PARS, QARS, QPRS, PIG32
Accession No.	Uniprot:P07814
Storage	6 months

## Application Details

Detection Range: 0.156-10ng/ml Assay Type: Sandwich Sensitivity: 0.055ng/ml

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.