

Human RC3H1 / Roquin-1 ELISA Kit

Catalog No: #EK21095



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human RC3H1 / Roquin-1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Roquin,RING finger and C3H zinc finger protein 1,RING finger and CCCH-type zinc finger domain-containing protein 1,RING finger protein 198,KIAA2025,RNF198
Accession No.	Uniprot:Q5TC82
Storage	6 months

Application Details

Detection Range:0.78-50ng/ml Assay Type:Sandwich Sensitivity:0.29ng/ml

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.