

# Rat Csnk1a1 / Casein kinase I isoform alpha ELISA Kit

Catalog No: #EK21096



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Csnk1a1 / Casein kinase I isoform alpha ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Other Names	CKI-alpha,CK1
Accession No.	Uniprot:P97633
Storage	6 months

## Application Details

Detection Range:78-5000pg/ml Assay Type:Sandwich Sensitivity:32pg/ml

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.