

Human DHRS3 / Short-chain dehydrogenase/reductase 3 ELISA Kit

Catalog No: #EK21097



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human DHRS3 / Short-chain dehydrogenase/reductase 3 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	DD83.1,Retinal short-chain dehydrogenase/reductase 1,retSDR1,Retinol dehydrogenase 17,Short chain dehydrogenase/reductase family 16C member 1,RDH17,SDR16C1,UNQ2424/PRO4983
Accession No.	Uniprot:O75911
Storage	6 months

Application Details

Detection Range:0.312-20ng/ml Assay Type:Sandwich Sensitivity:0.098ng/ml

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.