

# Human SSC5D / Soluble scavenger receptor cysteine-rich domain-containing protein SSC5D ELISA Kit



Catalog No: #EK21098

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human SSC5D / Soluble scavenger receptor cysteine-rich domain-containing protein SSC5D ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Soluble scavenger protein with 5 SRCR domains,SSc5D
Accession No.	Uniprot:A1L4H1
Storage	6 months

## Application Details

Detection Range:0.312-20ng/ml Assay Type:Sandwich Sensitivity:0.136ng/ml

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.