

Human BSG / Basigin ELISA Kit

Catalog No: #EK21099



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human BSG / Basigin ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	5F7, Collagenase stimulatory factor, Extracellular matrix metalloproteinase inducer, EMMPRIN, Leukocyte activation antigen M6, OK blood group antigen, Tumor cell-derived collagenase stimulatory factor, TCSF, CD147, UNQ6505/PRO21383
Accession No.	Uniprot:P35613
Storage	6 months

Application Details

Detection Range: 78-5000pg/mL Assay Type: Sandwich Sensitivity: 33pg/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.