

Human NPR1 / Atrial natriuretic peptide receptor 1 ELISA Kit

Catalog No: #EK21100



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human NPR1 / Atrial natriuretic peptide receptor 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Atrial natriuretic peptide receptor type A, ANP-A, Guanylate cyclase A, GC-A, ANPRA
Accession No.	Uniprot:P16066
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.15ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.