

# Human SCARB1 / Scavenger receptor class B member 1 ELISA Kit

Catalog No: #EK21101



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human SCARB1 / Scavenger receptor class B member 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	SRB1,CD36 and LIMPII analogous 1,CLA-1,CD36 antigen-like 1,Collagen type I receptor, thrombospondin receptor-like 1,SR-BI,CD36,CD36L1,CLA1
Accession No.	Uniprot:Q8WTV0
Storage	6 months

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.1ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.