

Human ACE2 / Angiotensin-converting enzyme 2 ELISA Kit

Catalog No: #EK21102



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ACE2 / Angiotensin-converting enzyme 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	ACE-related carboxypeptidase,Angiotensin-converting enzyme homolog,ACEH,Metalloprotease MPROT15,UNQ868/PRO1885
Accession No.	Uniprot:Q9BYF1
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.34ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.