

# Rat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit

Catalog No: #EK20211



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rat
Storage	12 months

## Application Details

Detect Range:15.6-1000pg/ml
Sensitivity:7.1pg/ml
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 &mu;L
Assay Time:1-4.5h
Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.