

## ZNF148 (Phospho-Ser306) Antibody

Catalog No: #SAB603

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ZNF148 (Phospho-Ser306) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified from rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	HM
Specificity	ACSS2 (phospho-Ser659) Antibody detects endogenous levels of ZNF148 only when phosphorylated at Ser306.
Immunogen Description	A synthesized peptide derived from human ZNF148 around the phosphorylation site of Ser306.
Other Names	Zinc finger protein 148; Transcription factor ZBP-89; Zinc finger DNA-binding protein 89; ZBP89
SDS-PAGE MW	89kDa
Concentration	1 mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000

## Background

Involved in transcriptional regulation. Represses the transcription of a number of genes including gastrin, stromelysin and enolase. Binds to the G-rich box in the enhancer region of these genes.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.