

Human Cytokeratin Fragment Antigen 21-1 (CYFRA21-1) ELISA Kit

Catalog No: #EK20139



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20139-1 48T #EK20139-2 96T

Description

Product Name	Human Cytokeratin Fragment Antigen 21-1 (CYFRA21-1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Hu

Application Details

Detect Range:0.312-20ng/ml

Sensitivity:0.124ng/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.