

Fish Thyrotropin Releasing Hormone (TRH) ELISAB Kit

Catalog No: #EK20148

Package Size: #EK20148-1 48T #EK20148-2 96T



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fish Thyrotropin Releasing Hormone (TRH) ELISAB Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Fish

Application Details

Detect Range: 0.5 - 100 pg/ml

Sensitivity: 0.4 pg/ml

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: competitive inhibition

Note: This product is for in vitro research use only and is not intended for use in humans or animals.