

# Recombinant 2019-nCoV S Protein RBD (N501Y, Mammalian, C-His)

Catalog No: #AG0069



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AG0069-1 100ug #AG0069-2 1mg

## Description

Product Name	Recombinant 2019-nCoV S Protein RBD (N501Y, Mammalian, C-His)
Purification	>90% as determined by SDS-PAGE
Other Names	RBD Protein, Spike RBD Protein, UK variant (B.1.1.7), Alpha
Calculated MW	35kDa
Formulation	Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

## Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.