

Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His)

Catalog No: #AG0071



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AG0071-1 100ug #AG0071-2 1mg

Description

| | |
|---------------|---|
| Product Name | Recombinant 2019-nCoV S Protein RBD (K417N, E484K, N501Y, Mammalian, C-His) |
| Purification | >90% as determined by SDS-PAGE |
| Other Names | RBD Protein, Spike RBD Protein, South African (ZA) variant (B.1.351), Beta |
| Calculated MW | 35kDa |
| Formulation | Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5 |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt. |

Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.