

Recombinant 2019-nCoV S Protein RBD (L452R, T478K, Mammalian, C-His)

Catalog No: #AG0073



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #AG0073-1 100ug #AG0073-2 1mg

Description

Product Name	Recombinant 2019-nCoV S Protein RBD (L452R, T478K, Mammalian, C-His)
Purification	>90% as determined by SDS-PAGE
Other Names	RBD Protein, Spike RBD Protein, IN variant (B.1.617.2), Delta
Calculated MW	35kDa
Formulation	Supplied as solution form in PBS, pH7.5 or lyophilized from PBS, pH7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week .Store at -20 to -80 °C for twelve months from the date of receipt.

Product Description

In general, proteins are shipped out with blue ice unless customers require otherwise.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.