

Sheep IgG

Catalog No: #IG301



Package Size: #IG301-1 10mg #IG301-2 1g

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sheep IgG
Purification	Affinity chromatography or caprylic acid-ammonium sulfate, β 95%
Concentration	β 10 mg/ml
Storage	-20°C

Product Description

Source: Normal, healthy, unimmunized sheep serum

Appearance: Liquid

Filtration: 0.2 μ m

Buffer :10mM sodium Phosphate150mM sodium chloridePBSPH7.4

Preservative: None

Stability: Stable for 24 months

Precautions: This product can be kept at 2-8 Λ C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.