

TNFAIP8L3 Polyclonal Antibody

Catalog No: #29583



Package Size: #29583-1 50ul #29583-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TNFAIP8L3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Immunogen Description	Synthesized peptide derived from human TNFAIP8L3.
Calculated MW	33kDa
Concentration	1mg/ml
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Images



Western Blot analysis of 1,293T 2, mouse-brain cells using primary antibody diluted at 1:1000(4°C overnight). Secondary antibody:Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour)

Background

caution:It is uncertain whether Met-1 or Met-89 is the initiator.,similarity:Belongs to the TNFAIP8 family.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.