

Anti-CD8/PE-Cy5.5 Conjugated Antibody

Catalog No: #C08440PE-Cy5.5



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anti-CD8/PE-Cy5.5 Conjugated Antibody
Brief Description	Anti-CD8/PE-Cy5.5 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Human,Mouse
Accession No.	Swiss_Prot:P01731 Gene_ID:12525
Calculated MW	27kDa
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow Cyt=3μg/Test

ICC=1:100-500

IF=1:100-500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.