

Recombinant Human DNA repair protein RAD50(RAD50) ,partial

Catalog No: #AP89629

Package Size: #AP89629-1 20ug #AP89629-2 100ug #AP89629-3 1mg



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DNA repair protein RAD50(RAD50) ,partial
Host Species	Yeast
Purification	>85% (SDS-PAGE)
Other Names	DNA repair protein RAD50; hRAD50; hsRAD50; Mrell; NBSLD; Rad 50; RAD50 2; rad50; RAD50 homolog (S cerevisiae); RAD50 homolog; RAD50 PEN; RAD50, S. cerevisiae, homolog of; RAD50_HUMAN; RAD502; Rad50I; Truncated RAD50 protein
Tag Info	N-terminal His-tagged
Target Sequence	1-235aa; partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.