

Recombinant Mouse Caveolin-1(Cav1)

Catalog No: #AP89632

Package Size: #AP89632-1 20ug #AP89632-2 100ug #AP89632-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Caveolin-1(Cav1)
Brief Description	Recombinant Mouse Caveolin-1(Cav1)
Target Name	His-tagged
Other Names	Cav1; Cav; Caveolin-1
Tag Info	N-terminal His-tagged
Target Sequence	2-178 SGGKYVDSE GHLYTVPIRE QGNIYKPNNK AMADEVTEKQ VYDAHTKEID LVNRDPKHLN DDVVKIDFED VIAEPEGTHS FDGIWKASFT TFTVTKYWFY RLLSTIFGIP MALIWGIYFA ILSFLHIWAV VPCIKSFLIE IQCISRVYSI YVHTFCDPLF EAIGKIFSNI RISTQKEI
Formulation	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.