

Human Angiostatin (ANG) ELISA Kit

Catalog No: #EK20212



Package Size: #EK20212-1 48T #EK20212-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human Angiostatin (ANG) ELISA Kit
Applications	ELISA
Species Reactivity	Hu

Application Details

Detect Range:78.125-5000pg/ml

Sensitivity:31pg/ml

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only and is not intended for use in humans or animals.