

# Ulex Europaeus Lectin 1 Polyclonal Antibody APC Conjugated

Catalog No: #C90390APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ulex Europaeus Lectin 1 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Ulex europaeus
Immunogen Type	101-200/243
Conjugates	APC
Target Name	Ulex Europaeus Lectin 1
Other Names	Ulex Europaeus Lectin 1
Accession No.	Uniprot:Q9BTC0Gene ID:
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.