

beta subunit Cholera Toxin Polyclonal Antibody Bio Conjugated

Catalog No: #C90889Bio



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	beta subunit Cholera Toxin Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Cholera Toxin
Conjugates	Bio
Target Name	beta subunit Cholera Toxin
Other Names	beta subunit Cholera Toxin
Accession No.	Uniprot:P01556Gene ID:
Calculated MW	12
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB □ 1:500-2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.