

CD80 Polyclonal Antibody PE Conjugated

Catalog No: #C91078PE



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD80 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC
Species Reactivity	Human, Mouse, Rat
Conjugates	PE
Target Name	CD80
Other Names	CD80
Accession No.	Uniprot:P33681Gene ID:941
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-CytoO 2µg/Test

Note: This product is for in vitro research use only and is not intended for use in humans or animals.