

Versican Polyclonal Antibody RBITC Conjugated

Catalog No: #C91752RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Versican Polyclonal Antibody RBITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IF |
| Species Reactivity | Ms Rat |
| Immunogen Type | 3301-3396/3396 |
| Conjugates | RBITC |
| Target Name | Versican |
| Other Names | Versican |
| Accession No. | Uniprot:P13611Gene ID:1462 |
| Calculated MW | 370 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.