

Ulex Europaeus Lectin 1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92007Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ulex Europaeus Lectin 1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Ulexeuropaeus
Immunogen Type	151-243/243
Conjugates	Bio
Target Name	Ulex Europaeus Lectin 1
Other Names	Ulex Europaeus Lectin 1
Accession No.	Uniprot:P22972Gene ID:
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.