

Vibrio Parahemolyticus Polar flagellin B/D Polyclonal Antibody Biotin Conjugated

Catalog No: #C92730Bio

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vibrio Parahemolyticus Polar flagellin B/D Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Vibrio Parahemolyticus
Immunogen Type	151-250/378
Conjugates	Bio
Target Name	Vibrio Parahemolyticus Polar flagellin B/D
Other Names	Vibrio Parahemolyticus Polar flagellin B/D
Accession No.	Uniprot:Q56702Gene ID:48370135
Calculated MW	40
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only and is not intended for use in humans or animals.