

Ferrous fumarate

Catalog No: #S0001

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ferrous fumarate
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	169.9
Formulation	C4H2FeO4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:Metabolism systemSMILES:C(=C/C(=O)[O-])\C(=O)[O-].[Fe+2]Pathways:OthersReceptor:OthersBoiling pt:Melting pt:>280CSolubility:Soluble in WaterAppearance:Brownness PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Garca-Casal MN, Leets I. J Food Sci. 2014 Apr;79(4):H706-12.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.