

Phenethyl alcohol

Catalog No: #S0011



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phenethyl alcohol
Brief Description	Inhibitors
Purification	99.82%
Target Name	Others
Calculated MW	122.16
Formulation	C8H10O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:ImmunologySMILES:OCCC1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:219-221CMelting pt:-27CSolubility:soluble in WaterAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lingappa BT, et al. Science. 1969 Jan 10;163(3863):192-4.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.