

Mecarbinat

Catalog No: #S0013

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mecarbinat
Brief Description	Inhibitors
Purification	99.12%
Target Name	Others
Calculated MW	233.27
Formulation	C13H15NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:OtherSMILES:CCOC(=O)C1=C(C)N(C)C2=C1C=C(O)C=C2Pathways:OthersReceptor:OthersBoiling pt:399.8CMelting pt:208-212CSolubility:DMSO: 47 mg/mL (201.49 mM)Appearance:Gray PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Shadurski? KS, et al. Farmakol Toksikol. 1983 Mar-Apr;46(2):115-20.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.