

cis-Anethol

Catalog No: #S0014

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	cis-Anethol
Brief Description	Inhibitors
Purification	98.29%
Target Name	Others
Calculated MW	148.2
Formulation	C10H12O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:Endocrine systemSMILES:C/C=C/C1=CC=C(C=C1)OCPathways:OthersReceptor:OthersBoiling pt:234-237CMelting pt:20-21CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Shimoni E, et al. J Biol Chem. 2002 Apr 5;277(14):11866-72.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.