

Mephenesin

Catalog No: #S0015

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mephenesin
Brief Description	Inhibitors
Purification	99.81%
Target Name	Others antagonist
Calculated MW	182.22
Formulation	C10H14O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:OtherSMILES:CC1=CC=CC=C1OCC(O)COPathways:OthersReceptor:OthersBoiling pt:153-154CMelting pt:69-71CSolubility:DMSO: 10 mMAppearance:White Cotton PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Prescire Int. 2013 May;22(138):127-8.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.