

Phenolphthalein

Catalog No: #S0017

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Phenolphthalein |
| Brief Description | Inhibitors |
| Purification | 99.09% |
| Target Name | Others |
| Calculated MW | 318.33 |
| Formulation | C20H14O4 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Product Description

Research Area: Metabolism

systemSMILES:OC1=CC=C(C=C1)C1(OC(=O)C2=CC=CC=C2)C1=CC=C(O)C=C1Pathways:OthersReceptor:OthersBoiling pt:548.7CMelting pt:258-263CSolubility:Water: 400 mg/LAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Uesawa Y, et al. Arch Toxicol. 2007 Mar;81(3):163-8.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.