

## Vanatone

Catalog No: #S0021

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vanatone
Brief Description	Inhibitors
Purification	99.38%
Target Name	Others
Calculated MW	150.17
Formulation	C9H10O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Product Description

Research Area:OtherSMILES:CCC(=O)C1=CC=C(O)C=C1Pathways:OthersReceptor:OthersBoiling pt:152-154CMelting pt:36-38CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Cizmrikov R, Kozlovsk? J. Ceska Slov Farm. 1994 Oct;43(5):226-8.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.