

Octisalate

Catalog No: #S0027

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Octisalate
Brief Description	Inhibitors
Purification	99.67%
Target Name	Others
Calculated MW	250.34
Formulation	C15H22O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:OtherSMILES:CCCC(CC)COC(=O)C1=CC=CC=C1OPathways:OthersReceptor:OthersBoiling pt:331.5CMelting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Singh M, Beck MH. Contact Dermatitis. 2007 Jan;56(1):48.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.