

Enoxolone

Catalog No: #S0036

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Enoxolone
Brief Description	Inhibitors
Purification	99.05%
Target Name	Dehydrogenase inhibitor
Calculated MW	470.68
Formulation	C30H46O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:Nervous

systemSMILES:C[C@]12CC[C@](C[C@H]1C3=CC(=O)[C@@H]4[C@]5(CC[C@@H](C([C@@H]5CC[C@]4([C@@]3(CC2)C)(C)C)O)(C)C(=O)O

Pathways:MetabolismReceptor:DHOBoling pt:588.3C?Melting pt:292-295CSolubility:DMSO: 94 mg/mL (199.71 mM)Appearance:White

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Boldt J, et al. Expert Opin PharmacOthers. 2007, 8(13):2135-47.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.