

Amidrine

Catalog No: #S0043



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Amidrine
Brief Description	Inhibitors
Purification	99.10%
Target Name	Others
Calculated MW	129.24
Formulation	C8H19N
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:Nervous systemSMILES:CC(C)CCCC(C)NPathways:OthersReceptor:OthersBoiling pt:155CMelting pt:Solubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. SEMENZA F. Boll Soc Ital Biol Sper. 1951 Mar;27(3 bis):354-7.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.