

## Corvotone

Catalog No: #S0048

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Corvotone
Brief Description	Inhibitors
Purification	99.86%
Target Name	GABA Receptor antagonist
Calculated MW	178.23
Formulation	C10H14N2O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Product Description

Research Area:Nervous systemSMILES:CCN(CC)C(=O)C1=CN=CC=C1Pathways:Membrane transporter/Ion channel;  
NeuroscienceReceptor:GABARBoiling pt:296-300CMelting pt:23CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility ,  
please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. ECKENHOFF JE, HAFKENSCHIEL JH. J Pharmacol Exp Ther. 1947 Dec;91(4):362-9.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.