

4-Aminoantipyrine

Catalog No: #S0050

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	4-Aminoantipyrine
Brief Description	Inhibitors
Purification	99.01%
Target Name	COX inhibitor
Calculated MW	203.24
Formulation	C11H13N3O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Product Description

Research Area:Nervous systemSMILES:CN1N(C(=O)C(N)=C1C)C1=CC=CC=C1Pathways:Immunology/Inflammation;

NeuroscienceReceptor:COXBoling pt:340 CMelting pt:105-110CSolubility:Soluble in WaterAppearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Brune K, Neubert A. Clin Exp Rheumatol. 2001 Nov-Dec;19(6 Suppl 25):S51-7.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.