

Oxeladin citrate

Catalog No: #S0056

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxeladin citrate
Brief Description	Inhibitors
Purification	99.12%
Target Name	Others
Calculated MW	527.6
Formulation	C26H41NO10
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Product Description

Research Area: Respiratory

systemSMILES:OC(=O)CC(O)(CC(O)=O)C(O)=O.CCN(CC)CCOCCOC(=O)C(CC)(CC)C1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:609.1C?Melting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Knopp C, Korsatko W. Pharmazie. 1982 Jan;37(1):76-7.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.