

## isoleucine

Catalog No: #S0063

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	isoleucine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Amino Acids and Derivatives inhibitor
Calculated MW	131.17
Formulation	C6H13NO2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Product Description

Research Area:Metabolism systemSMILES:CC[C@H](C)[C@@H](C(=O)O)NPathways:OthersReceptor:Amino Acids and DerivativesBoiling pt:Melting pt:288C?Solubility:soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Korman SH, et al. Clin Chem. 2005 Mar;51(3):610-7.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.